

# IV CANNULA WITH PERFORATED FLEXIBLE WINGS, WITH INJECTION PORT WITH SNAP FIT CAP

- Siliconized, Ultra sharp, Stainless Steel, back-cut needle to perform easy and painless venipuncture.
- Catheter available in FEP, PUR and PTFE material. Siliconised with proprietary technology ensuring smooth surface and minimal discomfort for both patient and clinician.
- Through the injection port it is possible to make injections at the same time as on-going fluid therapy. Snap fit cap for easy one-hand opening and closing.
- Large flashback Chamber allows quick visualization for successful venipuncture.
- Perforated wings provide ventilation & prevent skin maceration
- Transparent needle cover to prevent accidental damage to needle and catheter.

Size	Color code	Catheter O.D.mm	Catheter Length mm	Flow Rate (ml/mnt)	Box Qty (inner/outer)
14G	Orange	2.0	45	270	50/500
16G	Grey	1.7	45	200	50/500
17G	White	1.5	45	125	50/500
18G	Green	1.3	45	90	50/500
20G	Pink	1.1	32	60	50/500
22G	Blue	0.9	25	36	50/500
24G	Yellow	0.7	19	23	50/500



Find out more at [www.invumed.com](http://www.invumed.com)

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia  
invumed.com / info@invumed.com / +357 26220591  
© 2022 INVUMED Innovations LTD. All rights reserved. 16.12.22