

# PEDIATRIC IV CANNULA WITH SMALL WINGS, WITHOUT INJECTION PORT

- The Small wings & absence of injection port permits easy handling.
- Specially designed beveled back cut needle offers lower penetration force which is less painful to the patient.
- Siliconized, Ultra sharp, Stainless Steel, back-cut needle to perform easy and painless venipuncture.
- Catheter available in FEP, PUR and PTFE material. Siliconised with proprietary technology ensuring smooth surface and minimal discomfort for both patient and clinician.
- Large flashback chamber allows quick visualization for successful venipuncture.
- Transparent needle cover to prevent accidental damage to needle and catheter.

Size	Color code	Catheter O.D.mm	Catheter Length mm	Flow Rate (ml/mnt)	Box Qty (inner/outer)
24G	Yellow	0.7	19	23	50/500
26	Violet	0.6	19	17	50/500



Find out more at [www.invumed.com](http://www.invumed.com)

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia  
invumed.com / info@invumed.com / +357 26220591  
© 2022 INVUMED Innovations LTD. All rights reserved. 16.12.22

