

CVC KIT PREMIUM

CVC kit Premium has all the necessary sterile consumables for the CVC catheterisation in one ready-made kit provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- Introducer needle.
- J-tip guide wire with dispenser.
- 5 ml blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.
- Injection needle 22G.
- Injection needle 18G.
- Surgical Gloves.
- Surgical Pledgets.
- Surgery Sheet.
- Surgery Towel.
- Sterile Brushes.
- Gauze Pads.
- Suture with Needle.
- Wound Dressing.
- Scalpel.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB1P-14G20	14G	20 cm	0,81 mm × 60 cm	18G	8F
MB1P-16G20	16G	20 cm	0,81 mm × 60 cm	18G	6F
MB1P-18G20	18G	20 cm	0,50 mm × 60 cm	19G	5.5F
MB1P-20G10	20G	10 cm	0,40 mm × 50 cm	20G	4F
MB1P-20G13	20G	13 cm	0,40 mm × 50 cm	20G	4F
MB1P-20G15	20G	15 cm	0,40 mm × 50 cm	20G	4F
MB1P-22G10	22G	10 cm	0,20 mm × 50 cm	20G	4F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB2P-4F8	4F	8 cm	0,46 mm × 50 cm	20G	5.5F
MB2P-4F13	4F	13 cm	0,46 mm × 50 cm	20G	5.5F
MB2P-5F8	5F	8 cm	0,53 mm × 50 cm	20G	6F
MB2P-5F10	5F	10 cm	0,53 mm × 50 cm	20G	6F
MB2P-5F13	5F	13 cm	0,53 mm × 50 cm	20G	6F
MB2P-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB2P-8F20	8F	20 cm	0,89 mm × 60 cm	18G	9F
MB2P-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB3P-4.5F8	4.5F	8 cm	0,53 mm × 50 cm	20G	6F
MB3P-4.5F10	4.5F	10 cm	0,53 mm × 50 cm	20G	6F
MB3P-5.5F10	5.5F	10 cm	0,53 mm × 50 cm	20G	6.5F
MB3P-5.5F13	5.5F	13 cm	0,53 mm × 50 cm	20G	6.5F
MB3P-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB3P-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F

Quad lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB4P-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

