


CVC KIT OPTIMUM Y

CVC kit Optimum Y is sterile, pyrogen-free and atoxic and provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- Y-type introducer needle with check valve.
- J-tip guide wire with dispenser.
- 5 mL blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB10Y-14G15	14G	15 cm	0,81 mm × 50 cm	18G	8F
MB10Y-14G20	14G	20 cm	0,81 mm × 50 cm	18G	8F
MB10Y-14G30	14G	30 cm	0,81 mm × 70 cm	18G	8F
MB10Y-16G15	16G	15 cm	0,81 mm × 50 cm	18G	6F
MB10Y-16G20	16G	20 cm	0,81 mm × 50 cm	18G	6F
MB10Y-16G30	16G	30 cm	0,81 mm × 70 cm	18G	6F
MB10Y-18G15	18G	15 cm	0,50 mm × 50 cm	19G	5.5F
MB10Y-18G20	18G	20 cm	0,50 mm × 60 cm	19G	5.5F
MB10Y-18G30	18G	30 cm	0,50 mm × 70 cm	19G	5.5F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle Y-type	Dilator size
MB2OY-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB2OY-7F20	7F	20 cm	0,89 mm × 50 cm	18G	8F
MB2OY-8F15	8F	15 cm	0,89 mm × 50 cm	18G	9F
MB2OY-8F20	8F	20 cm	0,89 mm × 60 cm	18G	9F
MB2OY-8F30	8F	30 cm	0,89 mm × 70 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle Y-type	Dilator size
MB3OY-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB3OY-7F20	7F	20 cm	0,89 mm × 50 cm	18G	8F
MB3OY-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB3OY-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB3OY-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Quad lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle Y-type	Dilator size
MB4OY-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB4OY-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB4OY-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

