

INVUMED.COM

Infusion therapy

CENTRAL VENOUS CATHETER KITS

Invumed offers a comprehensive range of CVC kits for long-term central venous access for intensive infusion and/or transfusion therapy in different lengths, outside diameters, lumens and set configurations to meet both healthcare and patient demands.

- CVC kit Optimum
- CVC kit Optimum Y
- CVC kit Optimum ECG
- CVC kit Optimum Antimicrobial
- CVC kit Premium


CVC KIT OPTIMUM

CVC kit Optimum is sterile, pyrogen-free and atoxic.
Kit provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- Introducer needle.
- J-tip guide wire with dispenser.
- 5 ml syringe with injection needle.
- 5 ml blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB10-14G30	14G	30 cm	0,81 mm × 70 cm	18G	8F
MB10-16G30	16G	30 cm	0,81 mm × 70 cm	18G	6F
MB10-18G30	18G	30 cm	0,50 mm × 70 cm	19G	5.5F
MB10-20G30	20G	30 cm	0,40 mm × 70 cm	20G	4F
MB10-22G5	22G	5 cm	0,20 mm × 50 cm	20G	4F
MB10-22G8	22G	8 cm	0,20 mm × 50 cm	20G	4F
MB10-22G10	22G	10 cm	0,20 mm × 50 cm	20G	4F
MB10-22G13	22G	13 cm	0,20 mm × 50 cm	20G	4F
MB10-22G15	22G	15 cm	0,20 mm × 50 cm	20G	4F
MB10-22G16	22G	16 cm	0,20 mm × 50 cm	20G	4F
MB10-22G20	22G	20 cm	0,20 mm × 60 cm	20G	4F
MB10-22G30	22G	30 cm	0,20 mm × 70 cm	20G	4F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB20-4F5	4F	5 cm	0,46 mm × 50 cm	20G	5.5F
MB20-4F8	4F	8 cm	0,46 mm × 50 cm	20G	5.5F
MB20-4F10	4F	10 cm	0,46 mm × 50 cm	20G	5.5F
MB20-4F13	4F	13 cm	0,46 mm × 50 cm	20G	5.5F
MB20-4F15	4F	15 cm	0,46 mm × 50 cm	20G	5.5F
MB20-4F16	4F	16 cm	0,46 mm × 50 cm	20G	5.5F
MB20-4F20	4F	20 cm	0,46 mm × 60 cm	20G	5.5F
MB20-4F30	4F	30 cm	0,46 mm × 70 cm	20G	5.5F
MB20-5F30	5F	30 cm	0,53 mm × 70 cm	20G	6F
MB20-7F30	7F	30 cm	0,89 mm × 70 cm	18G	8F
MB20-8F30	8F	30 cm	0,89 mm × 70 cm	18G	9F
MB20-8.5F13	8.5F	13 cm	0,89 mm × 50 cm	18G	9F
MB20-8.5F16	8.5F	16 cm	0,89 mm × 50 cm	18G	9F
MB20-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB20-8.5F25	8.5F	25 cm	0,89 mm × 60 cm	18G	9F
MB20-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB30-4.5F5	4.5F	5 cm	0,53 mm × 50 cm	20G	6F
MB30-4.5F8	4.5F	8 cm	0,53 mm × 50 cm	20G	6F
MB30-4.5F10	4.5F	10 cm	0,53 mm × 50 cm	20G	6F
MB30-4.5F13	4.5F	13 cm	0,53 mm × 50 cm	20G	6F
MB30-4.5F15	4.5F	15 cm	0,53 mm × 50 cm	20G	6F
MB30-4.5F16	4.5F	16 cm	0,53 mm × 50 cm	20G	6F
MB30-4.5F20	4.5F	20 cm	0,53 mm × 60 cm	20G	6F
MB30-4.5F30	4.5F	30 cm	0,53 mm × 70 cm	20G	6F
MB30-5.5F30	5.5F	30 cm	0,53 mm × 70 cm	20G	6.5F
MB30-7F30	7F	30 cm	0,89 mm × 70 cm	18G	8F
MB30-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Quadruple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB40-8.5F5	8.5F	5 cm	0,89 mm × 50 cm	18G	9F
MB40-8.5F8	8.5F	8 cm	0,89 mm × 50 cm	18G	9F
MB40-8.5F10	8.5F	10 cm	0,89 mm × 50 cm	18G	9F
MB40-8.5F13	8.5F	13 cm	0,89 mm × 50 cm	18G	9F
MB40-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB40-8.5F16	8.5F	16 cm	0,89 mm × 50 cm	18G	9F
MB40-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB40-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22




CVC KIT OPTIMUM Y

CVC kit Optimum Y is sterile, pyrogen-free and atoxic and provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- Y-type introducer needle with check valve.
- J-tip guide wire with dispenser.
- 5 mL blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB10Y-14G15	14G	15 cm	0,81 mm × 50 cm	18G	8F
MB10Y-14G20	14G	20 cm	0,81 mm × 50 cm	18G	8F
MB10Y-14G30	14G	30 cm	0,81 mm × 70 cm	18G	8F
MB10Y-16G15	16G	15 cm	0,81 mm × 50 cm	18G	6F
MB10Y-16G20	16G	20 cm	0,81 mm × 50 cm	18G	6F
MB10Y-16G30	16G	30 cm	0,81 mm × 70 cm	18G	6F
MB10Y-18G15	18G	15 cm	0,50 mm × 50 cm	19G	5.5F
MB10Y-18G20	18G	20 cm	0,50 mm × 60 cm	19G	5.5F
MB10Y-18G30	18G	30 cm	0,50 mm × 70 cm	19G	5.5F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle Y-type	Dilator size
MB2OY-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB2OY-7F20	7F	20 cm	0,89 mm × 50 cm	18G	8F
MB2OY-8F15	8F	15 cm	0,89 mm × 50 cm	18G	9F
MB2OY-8F20	8F	20 cm	0,89 mm × 60 cm	18G	9F
MB2OY-8F30	8F	30 cm	0,89 mm × 70 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle Y-type	Dilator size
MB3OY-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB3OY-7F20	7F	20 cm	0,89 mm × 50 cm	18G	8F
MB3OY-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB3OY-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB3OY-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Quad lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle Y-type	Dilator size
MB4OY-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB4OY-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB4OY-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22




CVC KIT OPTIMUM ECG

CVC kit Optimum ECG is sterile, pyrogen-free and atoxic and provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- ECG cable allows for ECG monitoring of the catheter position by connecting to an adapter with subsequent signal conversion and transmission to an ECG monitor that displays the received data.
- Y-type introducer needle with check valve.
- J-tip guide wire with dispenser.
- 5 ml blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB10E-14G13	14G	13 cm	0,81 mm × 50 cm	18G	8F
MB10E-14G15	14G	15 cm	0,81 mm × 50 cm	18G	8F
MB10E-14G16	14G	16 cm	0,81 mm × 50 cm	18G	8F
MB10E-14G20	14G	20 cm	0,81 mm × 60 cm	18G	8F
MB10E-14G30	16G	30 cm	0,81 mm × 70 cm	18G	8F
MB10E-16G13	16G	13 cm	0,81 mm × 50 cm	18G	6F
MB10E-16G15	16G	15 cm	0,81 mm × 50 cm	18G	6F
MB10E-16G16	16G	16 cm	0,81 mm × 50 cm	18G	6F
MB10E-16G20	16G	20 cm	0,81 mm × 60 cm	18G	6F
MB10E-16G30	16G	30 cm	0,81 mm × 70 cm	18G	6F
MB10E-18G13	18G	13 cm	0,50 mm × 50 cm	19G	5,5F
MB10E-18G15	18G	15 cm	0,50 mm × 50 cm	19G	5,5F
MB10E-18G16	18G	16 cm	0,50 mm × 50 cm	19G	5,5F
MB10E-18G20	18G	20 cm	0,50 mm × 60 cm	19G	5,5F
MB10E-18G30	18G	30 cm	0,50 mm × 70 cm	19G	5,5F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB2OE-4F8	4F	8 cm	0,46 mm × 50 cm	20G	5.5F
MB2OE-4F10	4F	10 cm	0,46 mm × 50 cm	20G	5.5F
MB2OE-4F13	4F	13 cm	0,46 mm × 50 cm	20G	5.5F
MB2OE-4F15	4F	15 cm	0,46 mm × 50 cm	20G	5.5F
MB2OE-4F20	4F	20 cm	0,46 mm × 60 cm	20G	5.5F
MB2OE-5F8	5F	8 cm	0,53 mm × 50 cm	20G	6F
MB2OE-5F10	5F	10 cm	0,53 mm × 50 cm	20G	6F
MB2OE-5F13	5F	13 cm	0,53 mm × 50 cm	20G	6F
MB2OE-5F15	5F	15 cm	0,53 mm × 50 cm	20G	6F
MB2OE-5F20	5F	20 cm	0,53 mm × 60 cm	20G	6F
MB2OE-7F13	7F	13 cm	0,89 mm × 50 cm	18G	8F
MB2OE-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB2OE-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB2OE-7F30	7F	30 cm	0,89 mm × 70 cm	18G	8F
MB2OE-8F15	8F	15 cm	0,89 mm × 50 cm	18G	9F
MB2OE-8F20	8F	20 cm	0,89 mm × 60 cm	18G	9F
MB2OE-8F30	8F	30 cm	0,89 mm × 70 cm	18G	9F
MB2OE-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB2OE-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB2OE-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB3OE-4.5F8	4.5F	8 cm	0,53 mm × 50 cm	20G	6F
MB3OE-4.5F10	4.5F	10 cm	0,53 mm × 50 cm	20G	6F
MB3OE-4.5F13	4.5F	13 cm	0,53 mm × 50 cm	20G	6F
MB3OE-4.5F15	4.5F	15 cm	0,53 mm × 50 cm	20G	6F
MB3OE-4.5F20	4.5F	20 cm	0,53 mm × 60 cm	20G	6F
MB3OE-4.5F30	4.5F	30 cm	0,53 mm × 70 cm	20G	6F
MB3OE-5.5F8	5.5F	8 cm	0,53 mm × 50 cm	20G	6.5F
MB3OE-5.5F10	5.5F	10 cm	0,53 mm × 50 cm	20G	6.5F
MB3OE-5.5F13	5.5F	13 cm	0,53 mm × 50 cm	20G	6.5F
MB3OE-5.5F15	5.5F	15 cm	0,53 mm × 50 cm	20G	6.5F
MB3OE-5.5F20	5.5F	20 cm	0,53 mm × 60 cm	20G	6.5F
MB3OE-5.5F30	5.5F	30 cm	0,53 mm × 70 cm	20G	6.5F
MB3OE-7F13	7F	13 cm	0,89 mm × 50 cm	18G	8F
MB3OE-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB3OE-7F16	7F	16 cm	0,89 mm × 50 cm	18G	8F
MB3OE-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB3OE-7F30	7F	30 cm	0,89 mm × 70 cm	18G	8F

Quad lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB4AO-8.5F13	8.5F	13 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F16	8.5F	16 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB4AO-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22




CVC KIT OPTIMUM ANTIMICROBIAL

CVC kit Optimum Antimicrobial is sterile, pyrogen-free and atoxic and provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- Antimicrobial coating of Minocycline and Rifampicin effectively prevents the catheter related bloodstream infection (CRBSI).
- Introducer needle.
- J-tip guide wire with dispenser.
- 5 ml blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB1AO-14G15	14G	15 cm	0,81 mm × 50 cm	18G	8F
MB1AO-14G16	14G	16 cm	0,81 mm × 50 cm	18G	8F
MB1AO-14G20	14G	20 cm	0,81 mm × 50 cm	18G	8F
MB1AO-16G15	16G	15 cm	0,81 mm × 50 cm	18G	6F
MB1AO-16G16	16G	16 cm	0,81 mm × 50 cm	18G	6F
MB1AO-16G20	16G	20 cm	0,81 mm × 60 cm	18G	6F
MB1AO-16G30	16G	30 cm	0,81 mm × 70 cm	18G	6F
MB1AO-18G15	18G	15 cm	0,50 mm × 50 cm	19G	5.5F
MB1AO-18G16	18G	16 cm	0,50 mm × 50 cm	19G	5.5F
MB1AO-18G20	18G	20 cm	0,50 mm × 60 cm	19G	5.5F
MB1AO-18G30	18G	30 cm	0,50 mm × 70 cm	19G	5.5F
MB1AO-20G8	20G	8 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G10	20G	10 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G13	20G	13 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G15	20G	15 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G16	20G	16 cm	0,40 mm × 50 cm	20G	4F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB1AO-20G20	20G	20 cm	0,40 mm × 60 cm	20G	4F
MB1AO-20G30	20G	30 cm	0,40 mm × 70 cm	20G	4F
MB1AO-22G5	22G	5 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G8	22G	8 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G10	22G	10 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G13	22G	13 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G15	22G	15 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G16	22G	16 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G20	22G	20 cm	0,20 mm × 60 cm	20G	4F
MB1AO-22G30	22G	30 cm	0,20 mm × 70 cm	20G	4F

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB2AO-4F5	4F	5 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F8	4F	8 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F10	4F	10 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F13	4F	13 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F15	4F	15 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F16	4F	16 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F20	4F	20 cm	0,46 mm × 60 cm	20G	5.5F
MB2AO-5F8	5F	8 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F10	5F	10 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F13	5F	13 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F15	5F	15 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F16	5F	16 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F20	5F	20 cm	0,53 mm × 60 cm	20G	6F
MB2AO-5F30	5F	30 cm	0,53 mm × 70 cm	20G	6F
MB2AO-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB2AO-7F16	7F	16 cm	0,89 mm × 50 cm	18G	8F
MB2AO-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB2AO-7F30	7F	30 cm	0,89 mm × 70 cm	18G	8F
MB2AO-8F15	8F	15 cm	0,89 mm × 50 cm	18G	9F
MB2AO-8F16	8F	16 cm	0,89 mm × 50 cm	18G	9F
MB2AO-8F20	8F	20 cm	0,89 mm × 60 cm	18G	9F
MB2AO-8F30	8F	30 cm	0,89 mm × 70 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB3AO-4.5F8	4.5F	8 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F10	4.5F	10 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F13	4.5F	13 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F15	4.5F	15 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F16	4.5F	16 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F20	4.5F	20 cm	0,53 mm × 60 cm	20G	6F
MB3AO-4.5F30	4.5F	30 cm	0,53 mm × 70 cm	20G	6F
MB3AO-5.5F5	5.5F	5 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F8	5.5F	8 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F10	5.5F	10 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F13	5.5F	13 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F15	5.5F	15 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F16	5.5F	16 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F20	5.5F	20 cm	0,53 mm × 60 cm	20G	6.5F
MB3AO-5.5F30	5.5F	30 cm	0,53 mm × 70 cm	18G	6.5F
MB3AO-7F13	7F	13 cm	0,89 mm × 50 cm	18G	8F
MB3AO-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB3AO-7F16	7F	16 cm	0,89 mm × 50 cm	18G	8F
MB3AO-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB3AO-8.5F13	8.5F	13 cm	0,89 mm × 50 cm	18G	9F
MB3AO-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB3AO-8.5F16	8.5F	16 cm	0,89 mm × 50 cm	18G	9F
MB3AO-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB3AO-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Quad lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB4AO-8.5F13	8.5F	13 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F16	8.5F	16 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB4AO-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia

invumed.com / info@invumed.com / +357 26220591

© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22



CVC KIT PREMIUM

CVC kit Premium has all the necessary sterile consumables for the CVC catheterisation in one ready-made kit provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- Introducer needle.
- J-tip guide wire with dispenser.
- 5 ml blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.
- Injection needle 22G.
- Injection needle 18G.
- Surgical Gloves.
- Surgical Pledgets.
- Surgery Sheet.
- Surgery Towel.
- Sterile Brushes.
- Gauze Pads.
- Suture with Needle.
- Wound Dressing.
- Scalpel.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB1P-14G20	14G	20 cm	0,81 mm × 60 cm	18G	8F
MB1P-16G20	16G	20 cm	0,81 mm × 60 cm	18G	6F
MB1P-18G20	18G	20 cm	0,50 mm × 60 cm	19G	5.5F
MB1P-20G10	20G	10 cm	0,40 mm × 50 cm	20G	4F
MB1P-20G13	20G	13 cm	0,40 mm × 50 cm	20G	4F
MB1P-20G15	20G	15 cm	0,40 mm × 50 cm	20G	4F
MB1P-22G10	22G	10 cm	0,20 mm × 50 cm	20G	4F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB2P-4F8	4F	8 cm	0,46 mm × 50 cm	20G	5.5F
MB2P-4F13	4F	13 cm	0,46 mm × 50 cm	20G	5.5F
MB2P-5F8	5F	8 cm	0,53 mm × 50 cm	20G	6F
MB2P-5F10	5F	10 cm	0,53 mm × 50 cm	20G	6F
MB2P-5F13	5F	13 cm	0,53 mm × 50 cm	20G	6F
MB2P-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB2P-8F20	8F	20 cm	0,89 mm × 60 cm	18G	9F
MB2P-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB3P-4.5F8	4.5F	8 cm	0,53 mm × 50 cm	20G	6F
MB3P-4.5F10	4.5F	10 cm	0,53 mm × 50 cm	20G	6F
MB3P-5.5F10	5.5F	10 cm	0,53 mm × 50 cm	20G	6.5F
MB3P-5.5F13	5.5F	13 cm	0,53 mm × 50 cm	20G	6.5F
MB3P-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB3P-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F

Quad lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB4P-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

