


CVC KIT OPTIMUM ANTIMICROBIAL

CVC kit Optimum Antimicrobial is sterile, pyrogen-free and atoxic and provided with the following components:

- Single, double, triple or four lumen indwelling central venous catheter with soft tip, extension line, tube clamp and female Luer hubs.
- Antimicrobial coating of Minocycline and Rifampicin effectively prevents the catheter related bloodstream infection (CRBSI).
- Introducer needle.
- J-tip guide wire with dispenser.
- 5 ml blue wire introduction syringe.
- Vessel dilator 100 mm.
- Catheter holder and catheter holder clamp.

Single lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB1AO-14G15	14G	15 cm	0,81 mm × 50 cm	18G	8F
MB1AO-14G16	14G	16 cm	0,81 mm × 50 cm	18G	8F
MB1AO-14G20	14G	20 cm	0,81 mm × 50 cm	18G	8F
MB1AO-16G15	16G	15 cm	0,81 mm × 50 cm	18G	6F
MB1AO-16G16	16G	16 cm	0,81 mm × 50 cm	18G	6F
MB1AO-16G20	16G	20 cm	0,81 mm × 60 cm	18G	6F
MB1AO-16G30	16G	30 cm	0,81 mm × 70 cm	18G	6F
MB1AO-18G15	18G	15 cm	0,50 mm × 50 cm	19G	5.5F
MB1AO-18G16	18G	16 cm	0,50 mm × 50 cm	19G	5.5F
MB1AO-18G20	18G	20 cm	0,50 mm × 60 cm	19G	5.5F
MB1AO-18G30	18G	30 cm	0,50 mm × 70 cm	19G	5.5F
MB1AO-20G8	20G	8 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G10	20G	10 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G13	20G	13 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G15	20G	15 cm	0,40 mm × 50 cm	20G	4F
MB1AO-20G16	20G	16 cm	0,40 mm × 50 cm	20G	4F



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB1AO-20G20	20G	20 cm	0,40 mm × 60 cm	20G	4F
MB1AO-20G30	20G	30 cm	0,40 mm × 70 cm	20G	4F
MB1AO-22G5	22G	5 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G8	22G	8 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G10	22G	10 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G13	22G	13 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G15	22G	15 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G16	22G	16 cm	0,20 mm × 50 cm	20G	4F
MB1AO-22G20	22G	20 cm	0,20 mm × 60 cm	20G	4F
MB1AO-22G30	22G	30 cm	0,20 mm × 70 cm	20G	4F

Double lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB2AO-4F5	4F	5 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F8	4F	8 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F10	4F	10 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F13	4F	13 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F15	4F	15 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F16	4F	16 cm	0,46 mm × 50 cm	20G	5.5F
MB2AO-4F20	4F	20 cm	0,46 mm × 60 cm	20G	5.5F
MB2AO-5F8	5F	8 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F10	5F	10 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F13	5F	13 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F15	5F	15 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F16	5F	16 cm	0,53 mm × 50 cm	20G	6F
MB2AO-5F20	5F	20 cm	0,53 mm × 60 cm	20G	6F
MB2AO-5F30	5F	30 cm	0,53 mm × 70 cm	20G	6F
MB2AO-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB2AO-7F16	7F	16 cm	0,89 mm × 50 cm	18G	8F
MB2AO-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB2AO-7F30	7F	30 cm	0,89 mm × 70 cm	18G	8F
MB2AO-8F15	8F	15 cm	0,89 mm × 50 cm	18G	9F
MB2AO-8F16	8F	16 cm	0,89 mm × 50 cm	18G	9F
MB2AO-8F20	8F	20 cm	0,89 mm × 60 cm	18G	9F
MB2AO-8F30	8F	30 cm	0,89 mm × 70 cm	18G	9F

Triple lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB3AO-4.5F8	4.5F	8 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F10	4.5F	10 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F13	4.5F	13 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F15	4.5F	15 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F16	4.5F	16 cm	0,53 mm × 50 cm	20G	6F
MB3AO-4.5F20	4.5F	20 cm	0,53 mm × 60 cm	20G	6F
MB3AO-4.5F30	4.5F	30 cm	0,53 mm × 70 cm	20G	6F
MB3AO-5.5F5	5.5F	5 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F8	5.5F	8 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F10	5.5F	10 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F13	5.5F	13 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F15	5.5F	15 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F16	5.5F	16 cm	0,53 mm × 50 cm	20G	6.5F
MB3AO-5.5F20	5.5F	20 cm	0,53 mm × 60 cm	20G	6.5F
MB3AO-5.5F30	5.5F	30 cm	0,53 mm × 70 cm	18G	6.5F
MB3AO-7F13	7F	13 cm	0,89 mm × 50 cm	18G	8F
MB3AO-7F15	7F	15 cm	0,89 mm × 50 cm	18G	8F
MB3AO-7F16	7F	16 cm	0,89 mm × 50 cm	18G	8F
MB3AO-7F20	7F	20 cm	0,89 mm × 60 cm	18G	8F
MB3AO-8.5F13	8.5F	13 cm	0,89 mm × 50 cm	18G	9F
MB3AO-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB3AO-8.5F16	8.5F	16 cm	0,89 mm × 50 cm	18G	9F
MB3AO-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB3AO-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Quad lumen

REF	CVC catheter size	CVC catheter length	Guide wire	Introducer Needle	Dilator size
MB4AO-8.5F13	8.5F	13 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F15	8.5F	15 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F16	8.5F	16 cm	0,89 mm × 50 cm	18G	9F
MB4AO-8.5F20	8.5F	20 cm	0,89 mm × 60 cm	18G	9F
MB4AO-8.5F30	8.5F	30 cm	0,89 mm × 70 cm	18G	9F

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22

