

NEPHROSTOMY BALLOON DILATOR SETS

Nephrostomy balloon dilator is also known as a nephrotic balloon dilator. It's used to dilate the musculofascial tract and renal capsule establishing a nephrostomy tract during a percutaneous procedure.

- Designed to maintain 30 FR diameter under high pressure.
- Designed to provide effective, uniform dilation along the nephrostomy tract.
- Includes firm amplatz renal sheath.
- Includes Guide Wire: 0.035 inch.
- Rated for inflation pressure up to 17 ATM.
- Radiopaque marker to confirm balloon position.


Code key: NBD-XXYY-M

- NBD: Nephrostomy Balloon Dilator.
- XX: Balloon diameter, Fr.
- YY: Balloon length, cm.
- M: Inflation system (optional).
- S: 10cc syringe (optional).

Kit contains:

- High pressure nephrostomy balloon dilator.
- 30 Fr Amplatz Renal Sheath.
- High pressure 3 way stopcock.
- Guide Wire: 0.035 inch.

REF	Balloon size, Fr	Catheter size, Fr	Shaft length, cm	Syringe	Inflation device	Box Qty (inner/midbox/outer)
NBD-3012	30	7	55	—	—	1/15/100
NBD-3012-S	30	7	55	+	—	1/15/100
NBD-3012-M	30	7	55	—	+	1/15/100



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 16.01.23