



SUPRAPUBIC CYSTOSTOMY CATHETER KITS

Suprapubic catheterization refers to the placement of a drainage tube into the urinary bladder just above the pubic symphysis. This is typically performed for individuals who are unable to drain their bladder via the urethra.

Invumed offers wide range of suprapubic cystostomy catheter kits based on metallic or plastic splittable cannulas supplied with different types of polyurethane catheters or two ways foley catheters.

- The kits comes with metallic or plastic splittable cannula.
- The kit contains either: polyurethane single lumen catheter with J or malecot tip, with a metallic straightner inside, or 2 way latex or silicone balloon catheter.
- The kit contains all accessories that facilitate easy insertion and post operative closed system drainage.
- Single-use, supplied sterile in a single blister.







Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
invumed.com / info@invumed.com / +357 26220591
© 2022 INVUMED Innovations LTD. All rights reserved. 11.01.23



SUPRAPUBIC CYSTOSTOMY CATHETER KITS SIZE RANGE

Catheter size, Ch/Fr	Cannula Hub Color Code				
		J-tip catheter	Malecot catheter	Latex Foley catheter	Silicone Foley catheter
		REF	REF	REF	REF
Splitable cystostomy metallic cannula					
8 Ch/Fr	Blue	SUP-08-J-K	SUP-08-M-8-K	SUP-08-L-8-K	SUP-08-S-8-K
10 Ch/Fr	Grey	SUP-10-J-K	SUP-10-M-12-K	SUP-10-L-12-K	SUP-10-S-12-K
12 Ch/Fr	White	SUP-12-J-K	SUP-12-M-12-K	SUP-12-L-12-K	SUP-12-S-12-K
14 Ch/Fr	Green	SUP-14-J-K	SUP-14-M-12-K	SUP-14-L-12-K	SUP-14-S-12-K
Splitable cystostomy plastic cannula					
8 Ch/Fr	Light blue	CYST-08-J-K	CYST-08-M-K	CYST-08-L-K	CYST-08-S-K
10 Ch/Fr	Black	CYST-10-J-K	CYST-10-M-K	CYST-10-L-K	CYST-10-S-K
12 Ch/Fr	White	CYST-12-J-K	CYST-12-M-K	CYST-12-L-K	CYST-12-S-K
14 Ch/Fr	Green	CYST-14-J-K	CYST-14-M-K	CYST-14-L-K	CYST-14-S-K
16 Ch/Fr	Orange	CYST-16-J-K	CYST-16-M-K	CYST-16-L-K	CYST-16-S-K
18 Ch/Fr	Red	CYST-18-J-K	CYST-18-M-K	CYST-18-L-K	CYST-18-S-K
20 Ch/Fr	Yellow	CYST-20-J-K	CYST-20-M-K	CYST-20-L-K	CYST-20-S-K
22 Ch/Fr	Purple	CYST-22-J-K	CYST-22-M-K	CYST-22-L-K	CYST-22-S-K
24 Ch/Fr	Navy Blue	CYST-24-J-K	CYST-24-M-K	CYST-24-L-K	CYST-24-S-K
Catheter kit includes:		<ul style="list-style-type: none"> J-tip polyurethane catheter with straightner Splitable cystostomy cannula Connection Tube Scalpel#11 Fixation wings Stopcock Urine bag (optional) 	<ul style="list-style-type: none"> Malecot polyurethane catheter with straightner Splitable cystostomy cannula Connection Tube Scalpel#11 Fixation wings Stopcock Urine bag (optional) 	<ul style="list-style-type: none"> Two way latex foley catheter Splitable cystostomy cannula Scalpel#11 Urine bag.(optional) 	<ul style="list-style-type: none"> Two way silicone foley catheter Splitable cystostomy cannula Scalpel#11 Urine bag.(optional)

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 11.01.23