


# NEOFUSER

## MICROINFUSION PUMPS WITH SINGLE FLOW RATE

A single-use infusion pumps with predefined basal flow rate for continuous infusion. Light and easy for an unrestricted early mobility. Pumps with preset flow rate, suitable for application in regional anesthesia, wound infiltration, as well as in the post-operative pain therapy, among others.

- Pump volume: 60(S), 100(M), 200(I), 275(L), 450(X), 600(G) ml.
- Basal rates: 0.5, 1.0, 1.5, 2.0, 4.0, 5.0, 8, 10, 13, 20, 50, 100, 200, 250 ml/h.
- Not made with latex, BPA, or DEHP.
- With Luer-Lock patient connector with color coded flow rate control unit.
- Efficient particle and self-venting 0,2 µm filter.
- Pressure-Activated Valve for one-way drug delivery.
- Portable clip designed to easily mount on IV pole stands.
- The clamp is used to start and stop the flow.



Basal Flow Rate, ml/hr	Color code	Ref. Nominal Volume					
		60 ml	100 ml	200 ml	275 ml	450 ml	600 ml
0,5	■	S05	M05	I05	L05	X05	G05
1	■	S10	M10	I10	L10	X10	G10
1.5	■	S15	M15	—	L15	X15	G15
2	■	S20	M20	I20	L20	X20	G20
4	■	S40	M40	I40	L40	X40	G40
5	■	S50	M50	I50	L50	X50	G50
8	■	S80	M80	I80	L80	X80	G80
10	■	S100	M100	I100	L100	X100	G100
13	■	S130	M130	I130	L130	X130	G130
20	□	—	M200	I200	L200	X200	G200
50	■	—	M500	I500	L500	X500	G500
100	■	—	M1000	I1000	L1000	X1000	G1000
200	■	—	M2000	I2000	L2000	X2000	G2000
250	■	—	M2500	I2500	L2500	X2500	G2500

Find out more at [www.invumed.com](http://www.invumed.com)

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia  
 invumed.com / info@invumed.com / +357 26220591  
 © 2022 INVUMED Innovations LTD. All rights reserved. 22.12.22