


NEOFUSER PLUS

MICROINFUSION PUMPS WITH SINGLE FLOW RATE WITH PCA

Uniqueflow B+ is a Patient Controlled Analgesia (PCA) pump for the postoperative analgesia. With Uniqueflow B+ the pain medication is administered as a combination of the preset continuous basal rate and the patient controlled bolus administration.

- Pump volume: 60(S), 100(M), 200(I), 275(L), 450(X) ml.
- Basal rates: 0.5, 1.0, 2.0, 4.0, 5.0, 8, 10 ml/h.
- Bolus volumes: 0.5(R), 1.0(E), 2.0(W) ml.
- Bolus lock timer: 8, 15, 60 min.
- Not made with latex, BPA, or DEHP.
- With Luer-Lock patient connector.
- Efficient particle and self-venting 0,2 µm filter.
- Progressive flow indicator.
- Portable clip designed to easily mount on IV pole stands.
- The clamp is used to start and stop the flow of chemotherapy.



Basal Flow Rate, ml/hr	Z Bolus Volume, ml	YY Lock-out Period, min	Ref. Nominal Volume				
			60 ml	100 ml	200 ml	275 ml	450 ml
0,5	0.5(R), 1.0(E), 2.0(W)	08, 15, 60	S05YYZ	M05YYZ	I05YYZ	L05YYZ	X05YYZ
1	0.5(R), 1.0(E), 2.0(W)	08, 15, 60	S10YYZ	M10YYZ	I10YYZ	L10YYZ	X10YYZ
2	0.5(R), 1.0(E), 2.0(W)	08, 15, 60	S20YYZ	M20YYZ	I20YYZ	L20YYZ	X20YYZ
4	0.5(R), 1.0(E), 2.0(W)	08, 15, 60	S40YYZ	M40YYZ	I40YYZ	L40YYZ	X40YYZ
5	0.5(R), 1.0(E), 2.0(W)	08, 15, 60	S50YYZ	M50YYZ	I50YYZ	L50YYZ	X50YYZ
8	0.5(R), 1.0(E), 2.0(W)	08, 15, 60	S80YYZ	M80YYZ	I80YYZ	L80YYZ	X80YYZ
10	0.5(R), 1.0(E), 2.0(W)	08, 15, 60	—	—	—	—	X80YYZ

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 22.12.22