


NEOFUSER VARIO

MICROINFUSION PUMPS WITH MULTI-FLOW RATES

Neofuser Vario is single-use infusion pumps with multi-rate module offers the possibility to choose flow speed up to seven different rates. Light and easy for an unrestricted early mobility. Neofuser Vario is ideal for all fields in pain therapy including regional anaesthesia.

- Pump volume: 60(S), 100(M), 200(I), 275(L), 450(X), 600(G) ml.
- Flow rate: 0 – 2 – 4 – 6 – 8 – 10 – 12 – 14 ml/h.
- Not made with latex, BPA, or DEHP.
- With Luer-Lock patient connector.
- Efficient particle and self-venting 0,2 µm filter.
- Pressure-Activated Valve for one-way drug delivery.
- Portable clip designed to easily mount on IV pole stands.
- The clamp is used to start and stop the flow.



Basal Flow Rate, ml/hr	Ref. Nominal Volume					
	60 ml	100 ml	200 ml	275 ml	450 ml	600 ml
0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5	VS0535	VM05035	VI0535	VL0535	VX0535	—
1.0, 2.0, 3.0, 4.0, 5.0, 6.0, 7.0	VS17	VM17	VI17	VL17	VX17	VG17
2.0, 4.0, 6.0, 8.0, 10.0, 12.0, 14.0	VS214	VM214	VI214	VL214	VX214	VG214
1.0, 2.0, 3.0	VS13	VM13	VI13	VL13	VX13	—
2.0, 3.0, 5.0	VS25	VM25	VI25	VL25	VX25	—
2.0, 4.0, 6.0	VS26	VM26	VI26	VL26	VX26	—

Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 22.12.22