

# CIRCULAR STAPLERS FOR RECTAL PROLAPSE AND HAEMORRHOIDS

Circular Staplers for Rectal Prolapse and Haemorrhoids are used for surgical operation according to the Longo method.

- Two stepped rows of staples, staggered, create a closed suture in the form of a circle.
- A knife is integrated into the head for simultaneous tissue resection during stitching.
- Dual safety lock out mechanism prevents misfiring.
- Fixed anvil to eliminate risk of detachment from the cartridge.
- Incremental markings to clearly measure distance from the dentate line.
- Green firing indicator window.
- 4 suture threading ports for even selection of prolapsed mucosa.

Supplied with the following 5 accessories:

- Clear Circular Anal Dilator,
- Clear Butterfly Anal Dilator,
- Purse-String Suture Anoscope,
- Suture Threader.

## Circular Stapler for Rectal Prolapse and Haemorrhoids Size range

Color code of staples reload part	Working part diameter, mm	Staple height before stapling, mm	Stapling gap sizes, mm
WHITE	32	4.2	0.75-2.0
GREEN	32	4.2	0.75-2.0



Find out more at [www.invumed.com](http://www.invumed.com)

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia  
 invumed.com / info@invumed.com / +357 26220591  
 © 2022 INVUMED Innovations LTD. All rights reserved. 27.12.22