

HASSON TROCARS

Hasson trocars are designed for surgeons who prefer the security of the open technique for the primary port placement.

Trocars blunt tip design prevents inadvertent damage to viscera during insertion.

The cannula depth may be rapidly adjusted by means of the simple and effective fixation device - threaded Anchor or balloon.

Ref	Length mm	O.D. mm	Type of fixation	Package
1850-V10MHBR	100	10	Ring and balloon fixation	5
1850-V12MHBA	100	12	Threaded Anchor and balloon fixation	5
1850-V12MHBR	100	12	Ring and balloon fixation	5
1850-V11MH	100	11	Fixation Cannula	5
1850-V12MH	100	12	Fixation Cannula	5
1850-V15MH	100	15	Fixation Cannula	10



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia

invumed.com / info@invumed.com / +357 26220591

© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22