

# SMOKE FILTER

A passive smoke evacuation system removes particles, cells, virus and odor causing noxious chemicals from surgical smoke generated in laparoscopic procedures while maintaining a clear field of vision.

It is indicated for use during any minimally invasive surgery involving insufflation, electrocautery, laser, or ultrasonic scalpel use.

Three-layer non-woven membrane absorbs the smoke generated from coagulation without distorting visibility during the operation.

The universal Luer connector meets ISO standards for leak free connection to the trocar.

Ref	Size, mm	Package
1850-VFS	10	1 / 50 / 500



Find out more at [www.invumed.com](http://www.invumed.com)

INVUMED Innovations LTD. Cyprus, 23 John Kennedy Avenue «Globe House» - 1075 Nicosia  
invumed.com / info@invumed.com / +357 26220591  
© 2022 INVUMED Innovations LTD. All rights reserved. 23.12.22